

SonoTip[®] Pro Control

Endoscopic Ultrasound-Guided FNA System – Single Use



- **Twist-Lock Technology for Sheath & Needle Length Adjustment**
Single Hand Operation

Needle Length Adjustment (0 - 8.5 cm)

- **Ergonomically Shaped Sturdy Handle Design**
Effective and Easy Puncture in Dense Target Lesions

Optimized for a Wide Range of Puncture Techniques

- **Specially Treated Needle**
Superior Needle Visibility

- **Versions with Unique Sheath Stabilizer**
Controlled Needle Advancement in Large Channel Scopes

- **Colored Plastic Sheath**
Improved Visualization Under Endoscopic Imaging

More Precise Needle Placement and Better Performance in Angulated Endoscope Positions




 **Medi-Globe[®]**
The Spirit of Care.™

SonoTip® Pro Control


25 Gauge

SonoTip® Pro Control – Including a VacLock Aspiration Syringe with a Filling Volume of 20 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-33-18-025	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath	25 Gauge 0.5 mm	1.8 mm	
GUS-33-27-025	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm)	25 Gauge 0.5 mm	2.7 mm	

22 Gauge

SonoTip® Pro Control – Including a VacLock Aspiration Syringe with a Filling Volume of 20 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-33-18-022	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath → Compatible with Guide Wires of maximal Ø .018"	22 Gauge 0.7 mm	1.8 mm	
GUS-33-27-022	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm) → Compatible with Guide Wires of maximal Ø .018"	22 Gauge 0.7 mm	2.7 mm	

SonoTip® Pro Control – Including a VacLock Aspiration Syringe with a Filling Volume of 20 ml

Catalog No.	Description	Needle Diameter	Outer Diameter	Beveled Stylet
GUS-34-18-022	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath Compatible with Guide Wires of maximal Ø .018"	22 Gauge 0.7 mm	1.8 mm	
GUS-34-27-022	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath with Tip Stabilizer to Accommodate Large Working Channels (e.g. 3.2 mm, 3.7 mm, 3.8 mm) → Compatible with Guide Wires of maximal Ø .018"	22 Gauge 0.7 mm	2.7 mm	


19 Gauge

SonoTip® Pro Control – Including a VacLock Aspiration Syringe with a Filling Volume of 20 ml


Catalog No.	Description	Needle Diameter	Outer Diameter	Round Stylet
GUS-33-21-019	EUS-Guided FNA System with Twist-Lock Technology → Blue Colored Plastic Sheath → Compatible with Guide Wires of maximal Ø .030" (Recommended Medi-Globe catalog no. GGW-08-30-400)	19 Gauge 1.1 mm	2.1 mm	

Components for SonoTip® Pro Control EUS-FNA Line


Aspiration Syringe with a Larger Filling Volume of 20 ml

Catalog No.	Description	Aspiration Syringe
VAC120	VacLok Syringe 20 ml (included in SonoTip Pro Control® EUS FNA Systems)	

Nitinol Guide Wire Compatible to 19 Gauge EUS-FNA Systems


Catalog No.	Description	Diameter	Length	Nitinol Guide Wire
GGW-08-30-400	Cyst Guide Wire, uncoated with 15 cm Radiopaque Flexible Spring Tip	.030"	400 cm	

Pseudo-Cyst Drainage Stents, Double Pigtail with Position Markers

Catalog No.	Description	French Size	Length	Pseudo-Cyst Drainage Stent
GBS-11-08-002	Double Pigtail Stent with Position Markers	8.5	2 cm	
GBS-11-08-005		8.5	5 cm	
GBS-11-08-007		8.5	7 cm	
GBS-11-08-009		8.5	9 cm	
GBS-11-10-002		10.0	2 cm	
GBS-11-10-005		10.0	5 cm	
GBS-11-10-007		10.0	7 cm	
GBS-11-10-009		10.0	9 cm	

SonoTip[®] Pro Control

Endoscopic Ultrasound-Guided FNA System – Single Use



Distributed by:

Ver. 6.0 (2014-01)